

PCT

WORLD INTELLECTUAL PROPERTY ORGANIZATION
International Bureau



INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(51) International Patent Classification 6 : A01N 1/02
(11) International Publication Number: WO 97/16965
(43) International Publication Date: 15 May 1997 (15.05.97)

(21) International Application Number: PCT/GB95/02635
(22) International Filing Date: 9 November 1995 (09.11.95)
(81) Designated States: AU, CA, JP, US, European patent (AT, BE, CH, DE, DK, ES, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE).

Published
With international search report.

(71) Applicants (for all designated States except US): MEDICAL RESEARCH COUNCIL [GB/GB]; Centre for Reproductive Biology, Chalmers Street, Edinburgh EH3 9EW (GB).
APPLIED RESEARCH SYSTEMS ARS HOLDING N.V. [NL/NL]; 14 John B Gostinweg, Curaçao (AN).

(72) Inventor; and
(75) Inventor/Applicant (for US only): ATKEN, Robert, John [GB/GB]; Medical Research Council, Centre for Reproductive Biology, 37 Chalmers Street, Edinburgh EH3 9EW (GB).

(74) Agents: CHAPMAN, Paul, William et al.; Kilburn & Strode, 30 John Street, London WC1N 2DD (GB).

(54) Title: USE OF CYB MEDIUM FOR THE TRANSPORTATION AND STORAGE OF SPERM

(57) Abstract

The use of CYB medium for the storage and/or transportation of semen at ambient temperature.

BEST AVAILABLE COPY

FOR THE PURPOSES OF INFORMATION ONLY

Codes used to identify States party to the PCT on the front pages of pamphlets publishing international applications under the PCT:

AM	Armenia	GB	United Kingdom	MW	Malawi
AT	Austria	GE	Georgia	MX	Mexico
AU	Australia	GN	Ghana	NE	Niger
BB	Barbados	GR	Greece	NL	Netherlands
BE	Belgium	HU	Hungary	NO	Norway
BF	Burkina Faso	IE	Ireland	NZ	New Zealand
BG	Bulgaria	IT	Italy	PL	Poland
BJ	Benin	JP	Japan	PT	Portugal
BR	Brazil	KE	Kenya	RO	Romania
BV	Bolivia	KG	Kyrgyzstan	RU	Russian Federation
CA	Canada	KP	Democratic People's Republic of Korea	SD	Sudan
CF	Central African Republic	KR	Republic of Korea	SE	Sweden
CG	Congo	KZ	Kazakhstan	SG	Singapore
CH	Switzerland	LI	Liechtenstein	SI	Slovenia
CI	Cote d'Ivoire	LK	Sri Lanka	SK	Slovakia
CM	Cameroon	LR	Liberia	SN	Senegal
CN	China	LT	Lithuania	SS	Swaziland
CS	Czechoslovakia	LU	Luxembourg	TD	Chad
CZ	Czech Republic	LV	Latvia	TG	Togo
DE	Germany	MC	Monaco	TJ	Tajikistan
DK	Denmark	MD	Republic of Moldova	TT	Trinidad and Tobago
EE	Estonia	MG	Madagascar	UA	Ukraine
ES	Spain	ML	Mali	UG	Uganda
FI	Finland	MR	Mauritania	US	United States of America
FR	France			UZ	Uzbekistan
GA	Gabon			VN	Viet Nam

USE OF CYB MEDIUM FOR THE TRANSPORTATION AND STORAGE OF SPERM

The present invention relates to a novel method for storing and transporting sperm which ensures its viability.

Diagnosis of male infertility has become increasingly more complex over recent years in view of the increasing sophistication of andrological techniques. This has led to the creation of specialised andrology centres which possess the necessary expertise in techniques such as computerised image analysis, cell biology and biochemistry, and can thus offer the most up to date diagnostic techniques.

However, a limitation on the growth of such specialised centres is that patients have to attend the centre in order to provide a semen sample immediately prior to analysis. Fresh samples of semen are required since semen loses viability and functional competence when left in the presence of seminal plasma for any length of time, usually after 1 hour.

There thus exists a need to provide a method for storing/preserving semen samples at ambient temperature such that they can be sent to andrological centres for diagnostic testing, removing the need for patients to attend personally.

The present inventors have now found that semen can be stored at ambient temperature in a particular medium known as an egg yolk buffer.

Thus, in a first aspect, the present invention provides a medium for storage and/or transportation of semen at ambient temperature which comprises CYB medium.

CYB medium is a known cryoprotectant medium first described by Wendel, L. and Prins, G. S., J. Androl. 8: 41-47 (1987). It consists of the following components:

- 40% TES/TRIS buffer
- 30% sodium citrate/fructose solution
- 20% fresh egg yolk
- 1% pen-strep solution (10,000 IU).

TES = N'-Tris (hydroxymethyl)methyl-2-aminoethane sulphonic acid.

TRIS = Hydroxymethylaminomethane.

The TES/TRIS buffer contains 8.66 g TES and 2.06 g TRIS per 200 ml of water.

The sodium citrate/fructose solution contains 5.88 g sodium citrate and 4.0 g fructose per 200 ml of water.

The constituents of the CYB medium are mixed together and centrifuged at 600 g for 2 x 10 mins to remove particulate matter. The pH is adjusted to 7.4.

In a preferred embodiment, the medium also includes an antioxidant. Suitable antioxidants include α -tocopherol (vitamin E), catalase, glutathione and mannitol. It is believed that these components will act as free radical scavengers. A particularly preferred antioxidant is α -tocopherol, present at a concentration of 1mM.

In a second aspect, the invention provides CYB medium for use in the storage/transportation of semen at ambient temperature.

5 In a third aspect, the invention provides the use of CYB medium for the storage and/or transportation of semen at ambient temperature.

10 In a fourth aspect, the invention provides a method for storing and/or transporting semen at ambient temperature which comprises the step of bringing the semen into contact with CYB medium. In this aspect of the invention, the semen is diluted approximately 1:1 with the CYB medium.

15 In preferred embodiments of all of the above-noted aspects, the semen is human semen.

20 Preferred features of each aspect of the invention are as defined for each other aspect, *mutatis mutandis*.

The invention will now be described by way of the following examples which should not be construed as in any way limiting the invention.

EXAMPLE 1

(A) METHODS

30 (1) Sample Density Protocol

10 μ l of sperm suspension was added to 190 μ l of sperm diluting fluid (SDF; see below), mixed well, and both sides of the counting chamber were filled (improved

Neubauer ruling).

This was then left to settle.

5 The number of sperm in 5 large squares running diagonally across the slide was counted, each large square having 16 smaller squares.

10 When sperm heads lay on the border, only those that lay on the top and left side border were included.

The number of sperm in 5 large squares = million sperm per ml (10^6 /ml).

15 Sperm Diluting Fluid (SDF)

NaHCO₃ 50 g
Formalin 10 ml

Made up to 1 litre in distilled water.

20 (ii) Scoring Motility

25 10 μ l of sperm suspension is placed on a slide and covered with a 24x24mm coverslip. Using an eyepiece graticule with a square grid, an area is defined and motile sperm within that area are counted. Non-motile sperm within the area are counted. The field is moved and repeat process until at least 100 sperm have been counted. Motility (motile sperm as % of total sperm counted) is calculated.

30

(iii) Percoll Prepared Sperm

3 ml of 100% prepared percoll solution (see below) was placed in the bottom of a test tube. 3 ml of 50%

prepared percoll solution (see below) was carefully layered on top of the 100% percoll (50% percoll was layered 1ml at a time using a 1 ml pipette).

5 2 ml of semen sample was layered on top of the gradient column. This was then centrifuged at 1900 rpm for 20 mins (500 g).

10 Seminal plasma was removed from the top of the column (saved frozen).

15 The sperm was separated in bands (50% from the middle of the gradient at the 50/100 interface, and 100% from the bottom of the tube). This was resuspended in 7-10 ml of BWB (see below) and centrifuged at 1900 rpm for 5 mins.

The supernatant was removed and the sperm pellets resuspended in a known volume of BWB (500 ml- 1ml depending on recovery of sperm pellet).

20 Density was recorded and the sperm concentration was adjusted to $20 \times 10^6/\text{ml}$.

DISCONTINUOUS PERCOLL GRADIENTS

100% Percoll Solution - 100 ml

10 ml of 10 x Earle's Balanced Salts Solution
FLOW LABS, IRVINE, SCOTLAND

90 ml of percoll

PHARMACIA LKB BIOTECHNOLOGY AB, UPPSALA, SWEDEN

6 ml of albuminar 5%
ARMOUR PHARMACEUTICAL COMPANY, EASTBOURNE, ENGLAND

5 3 mg of sodium pyruvate
0.37 ml of sodium lactate
200 mg of sodium hydrogen carbonate (NaHCO_3)

50% Percoll Solution

10 100% percoll solution 1:1 diluted with BWB.

BWB PREPARATION

15 BWB Stock - Made up in 1 litre of distilled water.

20 NaCl 5.54 g
KCl 0.356 g
CaCl₂ (dihydrate) 0.250 g
KH₂PO₄ 0.162 g
MgSO₄·H₂O 0.294 g

25 BWB - 200 ml

30 NaHCO₃ 420 mg
Glucose 200 mg
Sodium pyruvate 6 mg
Albuminar 5% (= 0.3% Final) 12 ml
Sodium lactate 0.74 ml

35 Penicillin/streptomycin 2.0 ml
GIBCO

Hepes buffer 4.0 ml
FLOW LABORATORY, IRVINE, SCOTLAND

40 BWB stock 181 ml

45

(B) SAMPLE TREATMENT

(i) Patient semen samples were obtained in 30 ml sterile plastic sample containers. A 30 minute liquification period was allowed prior to recording sample volume, density and round cell count.

(ii) The samples were then split into two aliquots and treated as follows:

One aliquot was prepared on percoll gradients as detailed above.

The remaining aliquot was mixed with a fixed volume (4.0 ml) of CVB medium, thereby providing at least a 1:1 dilution for most semen samples. Sample motility was checked again before packaging in an insulated polystyrene box for despatch by a courier service over a period of 24 hours.

Aliquots prepared on Percoll gradients were analysed as follows:

(a) Acrosome Reaction Assay

Sperm samples were prepared according to the 50%/100% percoll gradient protocol, adjusting the sperm density to $20 \times 10^6/\text{ml}$.

200 μl of sperm was added to an equivalent volume of A23187 free acid to give a final concentration of 1.25 μM or 2.5 μM A23187.

The mixture was incubated at 37°C for 3 hours.

35

The mixture was washed at 500 g for 5 minutes and the supernatant removed and resuspended in BW at $20 \times 10^6/\text{ml}$.

Motility of the sample was recorded.

50 μl of sperm suspension was added to 500 μl of Hypo-osmotic Swelling Medium and incubated for 1 hour at 37°C.

The suspension was centrifuged for 5 minutes at 500 g, the supernatant discarded and the pellet resuspended in 50 μl of ice cold methanol, i.e. final sperm concentration $20 \times 10^6/\text{ml}$.

10 μl was put onto each spot of a Hendley slide and was allowed to air dry.

This was then overlaid with Lectin-Fluorescein Isothiocyanate (2 mg/ml in PBS) and incubated for 15 minutes in the dark.

Excess lectin was then washed with PBS.

25 A drop of Citifluor was put on each spot of the Hendley slide and the slide was then covered with a coverslip for examination under a fluorescence microscope.

Hypo-osmotic Swelling Medium

7.35 g sodium citrate
13.51 g fructose
1 litre of distilled water.

RESULTS

(a) Sperm motility

21 semen samples chosen from a random selection of donors were investigated. Motility before and after shipment in CYB medium was measured.

TABLE 1

Sperm Progressive Motility (MPO: a + b) in % motile on n = 21 donors		
Semen	t = 0	semen/CYB t = 24
8	11	24
11	17	14
17	22	16
22	26	20
26	30	46
30	32	65
32	36	30
36	44	30
44	47	48
47	48	59
48	49	55
49	53	58
53	54	30
54	60	39
60	62	56
62	64	71
64	73	68
73	45	64
45	76	55
76	80	61
80		69
mean = 46		mean = 47
SD = 22		SD = 19

These results show that Sperm Progressive Motility at t = 24 h in CYB is on average well related to the response at t = 0.

(b) Sperm Penetration Assay

Sperm select 1:1 was diluted with BMW.

The diluted sperm select was then loaded into flat capillary tubes (200 µm).

One end of the tube was capped with Critoseal (Mackay and Lynn, Edinburgh).

The open end of the sperm select capillary tubes was placed into the semen sample (50µl).

This was then incubated at 37°C for 30 min at an angle of 20° (approx).

The sperm select and Penetrak tubes were removed from the sperm, placed on a marked microscope tube and the number of sperm present at 1, 3 and 4.5 cm from the open end of the tube was counted (counting the number of sperm in 2 fields - using a x 40 objective).

Aliquots sent by courier

On arrival of the samples, the samples were analysed for motility loss before percoll preparation and diagnostic assays were repeated to determine if the transport diluent was capable of sustaining these important diagnostic parameters in the "field" situation. Experimental variation is reduced by the use of strict correlated protocols and by employing the same technician in the assessment of subjective assays such as the acrosome reaction test.

(b) Acrosome Reaction Assay

For the same 21 samples of semen as noted above, the sperm percentage of viable cells was as follows:

TABLE 2

Ionophore A23187	
Semen t = 3	semen/CYB t = 24
21	34
25	55
27	28
37	19
43	66
47	42
57	61
59	56
60	83
61	62
64	80
70	79
74	62
78	82
79	75
79	82
79	81
81	87
87	90
mean = 61 SD = 21	mean = 64 SD = 21

These data show that the % cells viable, following treatment with Ionophore A23187 at t = 24 h in CYB, is, on average, well related to the response at t = 3.

(c) Sperm Penetration Assay

Sperm penetration (at t = 1.5 cm, 3 cm and 4.5 cm) at t = 0 and after t = 24 h in CYB was as follows:

TABLE 3

Sperm penetration on n = 20 donors						
1.5cm t=0	3cm t=0	4.5cm t=0	1.5cm t=24	3cm t=24	4.5cm t=24	
251	34	2	166	27	8	
174	19	1	24	4	0	
121	6	0	35	7	0	
224	43	8	250	47	10	
110	19	2	85	18	1	
260	54	9	235	110	8	
196	20	3	358	74	6	
275	45	8	86	26	3	
170	31	25	284	99	28	
110	13	2	300	117	15	
2	0	0	3	1	0	
28	2	0	3	2	0	
192	45	9	186	43	1	
92	18	1	101	21	2	
114	11	0	109	12	3	
210	21	2	215	23	3	
300	40	7	270	38	8	
34	3	0	39	5	1	
170	37	5	156	32	6	
mean = 160 SD = 93	mean = 26 SD = 23	mean = 4 SD = 6	mean = 116 SD = 111	mean = 36 SD = 36	mean = 6 SD = 7	

There was no statistically significant difference between the values measured at t=0 and measured at t=24 hr in CYB.

Example 2

General Semen Analysis (other than motility) of Spermatozoa incubated in CYB for 24 hours

Semen liquefaction was assessed as being either non-existent, moderate or normal. In samples examined, the presence of mucous threads, ie incomplete liquefaction was also backed up by drawing the sample through a pipette to determine if the fluid flowed freely. After mixing 1:1 with CYB and incubation for 24 hours the assessment was repeated.

Sample Number	Semen t=0	Semen/CYB t=24hr
1	Normal	Normal
2	Normal	Normal
3	Moderate	Moderate
4	Normal	Normal
5	None	Moderate

5

Clearly, incubation in CYB has no significant effect on this observed property of sperm.

10

Abnormal form analysis was performed on sperm in semen and mixed 1:1 with CYB at t=0 hours and t=24 hours. The cells were fixed in standard formalin solution before analysis, the cells being observed using a x40 objective.

15

Sample Number	CYB t=0	CYB t=24	Semen t=0	Semen t=24	% diff t=0	% diff t=24
1	65	54	61	55	17	10
2	68	83	78	69	18	11
3	68	52	54	48	24	11
4	70	76	71	73	8	3
5	74	74	79	78	0	1

20

Clearly, incubation with CYB does not appear to alter sperm morphology. Any discrepancies can be assumed to be due to sampling error.

25

Measurement of Antisperm Antibodies in CYB

Samples were treated with the blood serum of patients previously known to have high antisperm antibody levels. Thus, antibody was transferred to the test samples, allowing a clear estimation of the effect of CYB on antibody detection after 24 hours.

5

The semen sample was allowed to liquefy for 30 minutes and was then prepared on mini percoll gradients by centrifugation at 600g for 5 minutes.

10

The semen was then washed twice in 0.4% BSA in Earles culture medium by centrifugation at 200g for 5 minutes. The semen was then resuspended in 500 μ l of 0.4% BSA in Earles culture medium.

15

500 μ l of donor semen was incubated with 500 μ l of serum for 1 hour. After incubation, the sample was washed twice with 0.4% BSA in Earles culture medium and then resuspended with 1ml of culture medium and then mixed 1:1 with CYB.

20

For the MAR test the washed sperm were reconstituted into their seminal plasma before the MAR scoring. This involved mixing the semen with O+ blood preparation and IgG anti-sera. The IBBA test is amore sensitive and specific assay for the detection of anti-sperm antibodies.

25

30

CLAIMS:

1. A medium for storage and/or transportation of semen at ambient temperature which comprises CYB medium.
2. CYB medium for use in the storage and/or transportation of semen at ambient temperature.
3. The use of CYB medium for the storage and/or transportation of semen at ambient temperature.
4. A method for the storage and/or transportation of semen at ambient temperature which comprises the step of bringing the semen into contact with CYB medium.
5. A method as claimed in claim 4 wherein the semen is diluted 1:1 with CYB medium.
6. A medium as claimed in claim 1 or claim 2, the use as claimed in claim 3, or the method as claimed in claim 4 or claim 5 wherein the CYB medium further comprises an antioxidant.
7. A medium, use or method as claimed in claim 6 wherein the antioxidant is α -tocopherol, catalase, glutathione or mannitol.
8. A medium, use or method as claimed in claim 7 wherein the antioxidant is α -tocopherol, present at a concentration of 1mM.
9. A medium, use or method as claimed in claim 7 or claim 8 wherein the semen is human semen.

	MAR t=0	MAR t=24 CYB	IBBA t=24 CYB
Control	Negative	Negative	IgA/IgG -ve
Stock IgG	100% positive	100% positive	IgA 50% +ve IgG 90% +ve
1 in 2 IgG	100% positive	100% positive	IgA 10% +ve IgG 90% +ve
1 in 10 IgG	100% positive	100% positive	IgA 10% +ve IgG 90% +ve

The IBBA results corresponded to the results gained by the serum samples in previous assays. Thus, the action of CYB does not appear to impair the detection of anti-sperm antibodies.

INTERNATIONAL SEARCH REPORT

International Application No
PL/GB 95/02635

A. CLASSIFICATION OF SUBJECT MATTER

IPC 6 AB1N1/02

According to International Patent Classification (IPC) or to both national classification and IPC

B. FIELDS SEARCHED

Minimum documentation searched (classification system followed by classification symbols)

IPC 6 AB1N

Documentation searched other than minimum documentation to the extent that such documents are included in the fields searched

Electronic data base consulted during the international search (name of data base and, where practical, search terms used)

C. DOCUMENTS CONSIDERED TO BE RELEVANT

Category	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
X	FERTILITY AND STERILITY, vol. 35, no. 2, February 1981, pages 285-288, XP000576874 D.G. JASKEY ET AL.: "Twenty-four to ninety-six-hour storage of human spermatozoa in TEST-yolk buffer." see abstract	1-9
X	FERTILITY AND STERILITY, vol. 34, no. 6, December 1980, pages 607-609, XP000576876 P.N. ZAVOS ET AL.: "Motility and enzyme activity of human spermatozoa stored for 24 hours at +5A C and -196A C." see page 607, column 2, line 26 - line 29 -/--	1-9

☒ Further documents are listed in the continuation of box C.

☐ Patent family members are listed in annex.

* Special categories of cited documents:

- 'A' document defining the general state of the art which is not considered to be of particular relevance
- 'E' earlier document published on or after the international filing date
- 'L' document which may throw doubt on priority claim(s) or which is cited to establish the publication date of another citation or other special reasons (as specified)
- 'O' document referring to an oral disclosure, use, exhibition or other means
- 'P' document published prior to the international filing date but later than the priority date claimed
- 'T' later document published after the international filing date which is cited to understand the principle or theory underlying the invention
- 'X' document of particular relevance; the claimed invention cannot be considered novel or cannot be considered to involve an inventive step when the document is taken alone
- 'Y' document of particular relevance; the claimed invention cannot be considered to involve an inventive step when the document is taken alone in view of the prior art
- 'N' document member of the same patent family

Date of the actual completion of the international search

1 August 1996

Date of mailing of the international search report

13.08.96

Name and mailing address of the ISA

European Patent Office, P.B. 5118 Patzlaan 2
NL-2200 HV Rijswijk
Tel. (+31-70) 340-2040, Tx. 31 651 epo nl,
Fax (+31-70) 340-3016

Authorized officer

Decorte, D

Form PCT/ISA/210 (second sheet) (July 1992)

INTERNATIONAL SEARCH REPORT

International Application No
PL/GB 95/02635

C. (Continuation) DOCUMENTS CONSIDERED TO BE RELEVANT

Category	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
X	JOURNAL OF DAIRY SCIENCE, vol. 72, no. 10, 1989, CHAPAIN, ILLINOIS pages 2683-2690, XP002089926 A. IJAZ ET AL.: "Induction of bovine sperm capacitation by TEST-yolk semen extender." see page 2684, column 1, paragraph 1	1-9
X	JOURNAL OF ANDROLOGY, vol. 8, no. 1, 1987, pages 41-47, XP000576950 L. WEIDEL ET AL.: "Cryosurvival of human spermatozoa frozen in eight different buffer systems." cited in the application see page 42, column 2, paragraph 1	1,2,9
A	CHEMICAL ABSTRACTS, vol. 94, no. 9, 1981 Columbus, Ohio, US; abstract no. 63088v, G. BRUECKNER ET AL.: "Effectors at storage of spermatozoa." XP002089927 see abstract & DERMATOL. MONATSSCHR., vol. 166, no. 12, 1988, pages 784-792	6-8

Form PCT/ISA/210 (continuation of second sheet) (July 1992)

**This Page is Inserted by IFW Indexing and Scanning
Operations and is not part of the Official Record**

BEST AVAILABLE IMAGES

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

☒ **BLACK BORDERS**

☐ **IMAGE CUT OFF AT TOP, BOTTOM OR SIDES**

☒ **FADED TEXT OR DRAWING**

☒ **BLURRED OR ILLEGIBLE TEXT OR DRAWING**

☐ **SKEWED/SLANTED IMAGES**

☐ **COLOR OR BLACK AND WHITE PHOTOGRAPHS**

☐ **GRAY SCALE DOCUMENTS**

☐ **LINES OR MARKS ON ORIGINAL DOCUMENT**

☒ **REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY**

☐ **OTHER:** _____

IMAGES ARE BEST AVAILABLE COPY.

As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.